

LED CHART PROJECTOR TCP-1000

HIGH CONTRAST PROJECTOR

DELIGHT IN SIGHT

Plug and project forever and
ever. No change of lamps.
No maintenance.



- Maintenance-free LED light
- Sophisticated design
- Programmable display sequence
- Brighter and clearer chart display
- Quick chart selection
- Cyrillic and DIN chart disc available




 **TOMEY**
TECHNOLOGY AND VISION

THE TOMEY LED CHART PROJECTOR TCP-1000



QUALITY IN DETAIL

The Tomey LED Chart Projector **TCP-1000** is designed for durability and efficient operation. Thanks to the cold LED light source the lamp is maintenance-free due to its up to 500.000 hours lifetime. The intuitive remote control allows you to program three customised display sequences. The **TCP-1000** will be delivered including remote control, projection screen and base for immediate use. You can also easily mount the system on your refraction unit or onto the wall.



SPECIFICATIONS

MEASUREMENT RANGES

| | |
|---------------------------------|--|
| Auto-off function | Adjustable 5 / 10 / 15 minutes |
| Working distance | 2.9 to 6.1 m |
| Projection magnification | 30 x |
| Projection tilt | approx. $\pm 15^\circ$ |
| Background illumination | 200 cd/m ² at 5 m (adjustable) |
| Binocular vision test | Red-green / polarising |
| Available charts | 36 |
| Masks | Single, horizontal, vertical, Red-green / polarising |
| Optotypes (standard) | Numbers, letters, pictures, Snellen E, Landolt |
| Optional disk | Cyrillic letters, DIN |

DIMENSIONS & ELECTRIC REQUIREMENTS

| | |
|--------------------------|------------------------|
| Dimensions WDH | 191 x 310 x 200 mm |
| Weight | 4 kg |
| Lamp | 3V LED (373 to 801 nm) |
| Input | AC 115 V to 230 V |
| Frequency | 50/60 Hz |
| Power consumption | 11 Watt |



Remote control

2011/10 - subject to change without notice

TOMEY EUROPE
TOMEY GmbH
Am Weichselgarten 19a
91058 Erlangen, Germany

Phone +49 9131 777 10, Fax +49 9131 777 1 20
Email info@tomey.de

TOMEY ASIA-PACIFIC
TOMEY CORPORATION JAPAN
2-11-33 Noritakeshinmachi
Nishi-ku, Nagoya 451-0051, Japan

Phone +81 52 581 5327, Fax +81 52 561 4735
Email intl@tomey.co.jp

TOMEY
TECHNOLOGY AND VISION

www.tomey.de